

Hausted® Converge Multi-Purpose Treatment and Exam Stretcher



Multi-task stretcher for
the ER, PACU, OR and
L&D Departments



Versatility

- **Stretcher and chair configurations** allow for multiple uses and help eliminate unneeded patient transfers
- **New Increased weight capacity of 720 lbs (evenly distributed)** designed to meet the needs of the bariatric patient
- **Integral, self-storing stirrups** allow quick and precise pelvic exams, labor assessment or emergency treatment
- **Telescoping IV pole included** (permanently mounted), additional IV/Accessory wells are at all four corners

Mobility

- **Positioning flexibility** from supine stretcher to a full chair - maximizing patient accessibility during examination
- **Patient Safety:** the built-in, self-storing footrest assists in appropriate patient positioning
- **Patient Comfort:** the leg section can be raised or lowered as needed
- **Brake-and-steer caster system** is accessible from all four corners
- **Two steer-lock casters** at foot to optimize better control

Product Specifications

MODELS	ELECTRIC/BATTERY POWERED 4E7-APR-ST 4D7-APR-ST*
Minimum Height	25.5 in [65 cm]
Maximum Height	33.5 in [85 cm]
Trendelenburg	12°
Patient Surface Dimensions	27 x 72 in [69 x 183 cm]
Overall Width	31.12 in [80 cm]
Overall Length	81 in [206 cm]
Side Rails	40 x 15 in [101 x 38 cm]
Fowler	0-80°
Leg Section	0-85°
Footrest	Self-storing
Bumper	4 corners
Casters	8 in Tente
Maximum Weight Capacity Evenly Distributed	720 lbs [327 kg]

