

C-arm Access of a Table.  
Mobility of a Stretcher.

Hausted® Fluoro-Track  
Fluoroscopy Capable  
Stretcher

*The Ideal Solution for Fluoroscopy,  
Pain Therapy, Endoscopy, Emergency,  
and Ambulatory Surgery.*



**HAUSTED**®  
PATIENT HANDLING SYSTEMS, LLC

**OPERATIONAL BENEFITS**

- The hourglass-shaped base with trimmed sides creates an excellent window for the C-arm
- The dual column base is easy to maneuver and provides superior stability for the patient
- The easy-to-operate sliding top enables a broader range of imaging capability without having to move the patient

**ECONOMIC BENEFITS**

- Improved case turnaround and streamlined patient handling because pre-indication, treatment, and recovery can be performed on the same stretcher
- Reduced potential for staff injury from patient transfer

**PATIENT OUTCOME AND SAFETY BENEFITS**

- Risk management and patient safety needs are addressed with a wider surface and higher rails than traditional fluoroscopy tables
- Improved quality of care and reduced risk of patient injury by eliminating multiple patient transfers



**GENERAL SPECIFICATIONS**

| MODEL                         | 493-RPC-ST                  |
|-------------------------------|-----------------------------|
| Minimum Height                | 23.5 in [60 cm]             |
| Maximum Height                | 36 in [91 cm]               |
| Trendelenburgs                | 15°                         |
| Patient Surface Dimensions    | 25.25 x 76 in [64 x 193 cm] |
| Overall Width                 | 29.75 in [76 cm]            |
| Overall Length                | 80.25 in [204 cm]           |
| Retracto Rails                | 15 x 64 in [38 x 162.5 cm]  |
| C-arm Clearance               | 23 in [58.5 cm]             |
| Casters                       | 8 in Tente [20 cm]          |
| Degree of Backrest Adjustment | 0-80°                       |
| Maximum Weight Capacity       | 500 lbs [227 kg]            |

*Sliding Top*  
12 in (30.5 cm) toward foot end  
6 in (15.2 cm) toward head end



← Slides 12 in.



← Centered →



Slides 6 in. →



PATIENT HANDLING SYSTEMS, LLC

334.215.5151 **Main**  
877.706.5151 **Toll Free**  
334.215.5150 **Fax**

2511 Midpark Drive  
Montgomery, AL 36109

[www.hausted.com](http://www.hausted.com)